

How to Use Your Respiclick®

- 1) Check the dose counter before each use.
- 2) Make sure the cap is closed before each dose is taken.
- 3) Hold the inhaler straight up. Pull the cap all the way back until you hear a “click”.
- 3) Stand or sit up straight. Breathe out all the way. DO NOT breathe into the inhaler.
- 4) Seal your lips tightly around the mouthpiece.
- 5) Breathe in FAST and DEEP (1-2 seconds) through your mouth. DO NOT block the vent above the mouthpiece.
- 6) When your lungs are full, take the inhaler out of your mouth.
- 7) Hold your breath for 10 seconds or as long as you can.
- 8) Breathe out normally.
- 9) Check the dose counter to make sure you breathed in the dose.
- 10) Close the cap over the mouthpiece after each use.
- 11) If you need another dose, repeat steps 1-10.



*The Respiclick® is a dry powder inhaler.
When you breathe in a dose,
you may or may not feel or taste
the medicine.*

Respiclick® Tips

Do I need to prime the Respiclick®?

- This inhaler does not need to be primed.

How do I care for my Respiclick®?

- Store this inhaler in the foil pouch and only open when you are ready to use it.
- Keep this inhaler in a dry place.
- Keep this inhaler at room temperature.
- Do not wash or take this inhaler apart.
- If the mouthpiece needs to be cleaned, wipe off with a dry cloth.

How do I know when I need a refill?

- The dose counter shows the number of doses left in the inhaler.
- When this inhaler is brand new, it has 200 doses.
- The dose counter counts down by 2.
- When the dose counter reaches 20, it will turn red. Contact the pharmacy to get a refill.
- When the dose counter reaches “0”, it is empty. Throw it away.

**Talk to your health care provider
if you have questions
about your Respiclick®.**